

Product Datasheet

RAD9B antibody (orb417179)

Description	RAD9B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Cell cycle checkpoint control protein RAD9B protein (274-408AA)
Target	RAD9B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6WBX8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.