


Product Datasheet

USP37 antibody (orb417239)

| | |
|----------------------------|--|
| Description | USP37 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC |
| Immunogen | Recombinant Human Ubiquitin carboxyl-terminal hydrolase 37 protein (671-797AA) |
| Target | USP37 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q86T82 |
| Dilution Range | IHC-P:1:200-500 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemistry of paraffin-embedde...