


Product Datasheet

HIST1H3A (phospho-Y41) antibody (orb417598)

| | |
|----------------------------|---|
| Description | HIST1H3A (phospho-Y41) antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF |
| Immunogen | Peptide sequence around site of Phospho-Tyr (41) derived from Human Histone H3.1 |
| Target | HIST1H3A |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P68431 |
| Dilution Range | IF:1:50-200 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence staining of HeLa cell...