

Product Datasheet

MFAP4 antibody (FITC) (orb417825)

| | |
|------------------------|---|
| Description | MFAP4 antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | FITC |
| Immunogen | Recombinant Human Microfibril-associated glycoprotein 4 protein (100-191AA) |
| Target | MFAP4 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P55083 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |