


Product Datasheet

SOX8 antibody (orb417951)

Description	SOX8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Transcription factor SOX-8 protein (292-386AA)
Target	SOX8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P57073
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of U87 whole cell...