

Product Datasheet

BZLF2 antibody (FITC) (orb418065)

Description	BZLF2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Epstein-Barr virus Glycoprotein 42 protein (34-223AA)
Target	BZLF2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03205
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.