

Product Datasheet

Histone H2AX (Phospho-S139) antibody (orb418647)

Description	Rabbit monoclonal antibody to Histone H2AX (Phospho-S139)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthesized peptide
Target	H2AFX
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P16104
Dilution Range	Recommended dilution: WB:1:3000-1:10000
Expiration Date	12 months from date of receipt.