

Product Datasheet

Acetyl-Histone H2B type 1-B (K20) antibody (orb418671)

Description	Rabbit monoclonal antibody to Acetyl-Histone H2B type 1-B (K20)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	A synthesized peptide
Target	HIST1H2BB
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P33778
Dilution Range	WB:1:500-1:2000, IF/ICC:1:50-1:500
Expiration Date	12 months from date of receipt.