

## Product Datasheet

### Human CRBN protein (orb419252)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Recombinant Human cereblon  |
| <b>Endotoxins</b>        | Not test.   |
| <b>Preservatives</b>     | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.  |
| <b>Form/Appearance</b>   | Liquid or Lyophilized powder  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Tag</b>               | N-terminal 6xHis-tagged   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Tag Info: N-terminal 6xHis-tagged<br>Expression Region: 1-442aa<br>Sequence Info: Full Length   |
| <b>Protein Sequence</b>  | MAGEGDQQDAAHNMGNHLPLLPAAESEEEDEMEVEDQDSKEAKKPNIIINFDTSL<br>PTSHTYLGADMEEFHGRTLHDDDDSCQVIPVLPQVMMILIPGQTLPLQLFHPQEV<br>MVRNLIQKDRTFAYLAYSINVQEREAQFGTTAEIYAYREEQDFGIEIVKVKAIQRQ<br>FKVLELRTQSDGIQAKVQILPECVLPSTMSAVQLESLNKCQIFPSKPVSREDQCS<br>YKWWQKYQKRKFHCANLTSWPRWLYSLYDAETLMDRIKKQLREWENLKDDSD<br>LPSNPIDFSYRVAACLPIDDVLRQLLKGSAIQLRCELDMNKCTSLCCKQCQET<br>EITTKNEIFSLSLCGPMAAYVNPVHGYVHETLTVYKACNLNLIGRPSTESHVFPGYA<br>WTVAQCKICASHIGWKFTATKKDMSPQKFWGLTRSALLPTIPDTEDEISPKVILC<br>L |
| <b>Purity</b>            | Greater than 90% as determined by SDS-PAGE.   |
| <b>Source</b>            | in vitro E.coli expression system   |
| <b>MW</b>                | 54.5 kDa  |
| <b>Uniprot ID</b>        | <a href="#">Q96SW2</a>  |
| <b>Expiration Date</b>   | 6 months from date of receipt.  |