

Product Datasheet

Human OPTN protein (orb419284)

Description	Recombinant Human Optineurin
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 9xHis-tagged
Note	For research use only
Application notes	Tag Info: C-terminal 9xHis-tagged Expression Region: 1-577aa Sequence Info: Full Length of Mature Protein
Protein Sequence	MSHQPLSCLTEKEDSPSESTGNGPPHLAHPNLDFTPEELLQMKELLTENHQL KEAMKLNQAMKGRFEELSAWTEKQKEERQFFEIQSKEAKERLMALSHENEKLL EELGKLGKSERSSDPTDDSRLPRAEAEQEKDQLRTQVVRLQAEKADLLGIVSE LQLKLNSSGSSSEDSFVEIRMAEGEAEGSVKEIKHSPGPTRTVSTGTALSKYRSRS ADGAKNYFEHEELTVSQQLLCLREGNQKVERLEVALKEAKERVSDFEKKTNSRSE IETQTEGSTEKENDEEKGPETVGVSEVEALNLQVTSFLKELQEAHTKLSEAELMKK RLQEKQALERKNSAIPSELNEKQELVYTNKKLELQVESMLSEIKMEQAKTEDEK SKLTVLQMTHNKLLQEHNNALKTIEELTRKESEKVDRAVLKELSEKLELAEKALAS KQLQMDQEMKQTIKQEEDELETMTILRAQMEVYCSDFHAERAAREKIHEEKEQLA LQLAVLLKENDAFEDGGRQSLMEMQSRHGARTSDSDQAYLVQRGAEDRDW RQQRNIPHIHSCPCKCGEVLDPDIDTLQIHVMDCII
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Baculovirus
MW	67.9 kDa
Uniprot ID	Q96CV9
Expiration Date	6 months from date of receipt.