


Product Datasheet

Mouse IL36G protein (orb419313)

Description	Recombinant Mouse Interleukin-36 gamma active
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 20 mM Tris-HCl, pH 8.5, 150mM NaCl with Tween-20
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-tagged Expression Region: 13-164aa Sequence Info: Full Length of Mature Protein
Protein Sequence	GRETPDFGGEVFDLDDQQVWIFRNQALVTVPRSHRVTPVSVTILPCKYPESLEQDKG IAIYLGIQNPDKCLFCKEVNGHPTLLKKEEKILDLYHHPEPMKPFLFYHTRTGGTST FESVAFPGHYIASSKTGNPIFLTSSKKGYYNINFNLDIKS
Purity	> 98% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	17.3 kDa
Uniprot ID	Q8R460
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse IL36G protein

