


Product Datasheet

Rat XCL1 protein (orb419316)

Description	Recombinant Rat Lymphotoctin active
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered PBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-tagged Expression Region: 22-114aa Sequence Info: Full Length of Mature Protein
Protein Sequence	VGTEVLQESICVSLRTQRLPVQKIKTYTIKEGAMRAVIFVTKRGLRICADPQAKWV KTAIKTV DGRASASKSKAETIPTQAQRSASTAVTLTG
Purity	> 97% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	10.0 kDa
Uniprot ID	P51672
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat XCL1 protein

