


Product Datasheet

Human SPRC protein (orb419320)

Description	Recombinant Human SPARC active
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered PBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-tagged Expression Region: 18-297aa Sequence Info: Full Length of Isoform PIGF-2
Protein Sequence	APQQEALPDETEVVEETVAEVTESVSGANPVQVEVGEFDDGAEETEEVVVAENP CQNHHCXGKVCELDENNTPMCVCQDPTSCPAPIGEFEKVCSDNKTDFDSSCH FFATKCTLEGTKKGHKLHLDYIGPCKYIPCLDSELTFFPLMRDWLKNVLVTLYE RDEDNLLTEKQKLRVKKIHENEKRLEAGDHPVELLARDFEKNYNMYIFVHWQ FGQLDQHPIDGYLSHTELAPLRAPLIPMEHCTTRFFETCDLDNDKYIALDEWAGC FGIKQKDIDKDLVI
Purity	> 98% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	32.7 kDa
Uniprot ID	P09486
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human SPRC protein