

Product Datasheet

Human Myoglobin protein (orb419965)

Description	Recombinant Human Myoglobin is produced by our E.coli expression system and the target gene encoding...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl,200mM NaCl,1mM DTT,20% Glycerol,pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	C-6His
Uniprot ID	P02144
Expiration Date	6 months from date of receipt.