

Product Datasheet

Monkey CD38 protein (orb420018)

Description	Recombinant Cynomolgus ADP-ribosyl Cyclase/cyclic ADP-ribose Hydrolase 1 is produced by our Mammalia...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	C-6His
Uniprot ID	Q5VAN0
Expiration Date	6 months from date of receipt.