

Product Datasheet

5-HT2B Receptor antibody (orb420038)

Description

5-HT2B Receptor antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human 5-HT2B. Please contact us for the exact immunogen sequence.

Target

5-HT2B Receptor

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

54 kDa

Uniprot ID

[Q02152](#), [P41595](#), [P30994](#)

NCBI

[000867](#), [000858](#)

Entrez

[3357](#)

Dilution Range


WB: 1:200-1000

Expiration Date

12 months from date of receipt.



Western blot analysis of rat adrenal gla...



IHC-P staining of human breast cancer ti...