


## Product Datasheet

### Human Prolactin R / PRLR Protein (orb420092)

<b>Description</b>	Human Prolactin R, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 25 - Asp 2...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Prolactin R
<b>Preservatives</b>	50 mM Tris, 150 mM NaCl, pH7.5
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.3 kDa. The protein migrates as 32-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Purity</b>	95%
<b>Source</b>	Human Prolactin R, His Tag (orb420092) is expressed from human 293 cells (HEK293). It contains AA Gln 25 - Asp 234 (Accession # P16471-1).
<b>MW</b>	26.3 kDa
<b>Uniprot ID</b>	<a href="#">P16471</a>
<b>NCBI</b>	<a href="#">NP_000940.1</a>
<b>Expiration Date</b>	6 months from date of receipt.



Immobilized  
 Human  
 Prolactin,  
 Mouse  
 IgG2a...