

Product Datasheet

Mouse VEGF R1 / Flt-1 Protein (orb420098)

Description

Mouse VEGF R1 Protein, Mouse IgG2a Fc Tag is

AA...	AA...
Reactivity	Mouse
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	VEGF R1
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-mFc(mlgG2a)
Note	For research use only
Application notes	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 109.8 kDa. The protein migrates as 140 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	85%
Source	Mouse VEGF R1 Protein, Mouse IgG2a Fc Tag (orb420098) is expressed from human 293 cells (HEK293). It contains AA Tyr 23 - Glu 759 (Accession # P35969-1).
MW	109.8 kDa
Uniprot ID	P35969
NCBI	NP_034358.2
Expiration Date	6 months from date of receipt.



Loaded
Biotinylated
Mouse
VEGF164,
His, ...