
Product Datasheet

**F(ab')₂ Rat IgG F(ab')₂ antibody (Peroxidase)
(orb420167)**

Description

 F(ab')₂ Rat IgG F(ab')₂ antibody

Species/Host	Goat
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Immunogen	Anti-Rat IgG F(ab') ₂ fragment was produced by repeated immunization with Rat IgG F(ab') ₂ fragment in goat.
Preservatives	None
Form/Appearance	Lyophilized
Concentration	0.8 mg/mL
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Suitable for immunoblotting (western or dot blot), ELISA, immunoperoxidase electron microscopy and immunohistochemistry as well as other peroxidase-antibody based enzymatic assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Isotype	IgG F(ab') ₂
Clonality	Polyclonal
Purity	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

 reaction was observed against anti-
Pepsin, anti-Goat IgG F(c), Rat IgG F(c)