

Product Datasheet

TOPK (Phospho-Thr9) antibody (orb420854)

Description	Rabbit polyclonal antibody to TOPK (Phospho-Thr9)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	PBK
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36085
Uniprot ID	Q96KB5
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.