

Product Datasheet

eEF2k (Phospho-Ser366) antibody (orb420859)

Description	Rabbit polyclonal antibody to eEF2k (Phospho-Ser366)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	EEF2K
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	82144
Uniprot ID	O00418
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.