

Product Datasheet

Raf1 (Phospho-Ser296) antibody (orb420922)

Description	Rabbit polyclonal antibody to Raf1 (Phospho-Ser296)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	RAF1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	73052
Uniprot ID	P04049
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.