

Product Datasheet

Nipah Virus G Protein Antibody (orb421078)

Description	Monoclonal Antibody to Nipah G Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant ectodomain of Nipah G Protein
Target	No cross reaction with Dengue, Zika and Chikungunya envelope proteins
Preservatives	No Preservative
Form/Appearance	Sterile Filtered Phosphate Buffered Saline, pH 7.4
Concentration	Batch dependant
Storage	Store at 2-8°C
Note	For research use only
Application notes	Pairing for Immunoassay: Capture: orb421078 - Detection: orb421077
Isotype	IgG1
Clonality	Monoclonal
Purity	> 85% pure (SDS-PAGE). Protein A Chromatography
Dilution Range	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application
Expiration Date	12 months from date of receipt.