



## Product Datasheet

# Cardiac FABP Fabp3 Antibody (monoclonal, 6C4) (orb421121)


<b>Description</b>	Cardiac FABP Fabp3 Antibody (monoclonal, 6C4)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	E. coli-derived rat Cardiac FABP recombinant protein (Position: A2-A133).
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Storage</b>	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	15 kDa
<b>Uniprot ID</b>	<b>P07483</b>
<b>Expiration Date</b>	12 months from date of receipt.



WB analysis of Cardiac FABP using anti-C...



Sandwich ELISA - Recombinant rat Cardiac...



IHC analysis of Cardiac FABP using anti-...