

## Product Datasheet

### SR9011 (orb422178)

<b>Description</b>	SR9011 is a REV-ERB $\alpha$ / $\beta$ agonist with IC50s of 790 nM and 560 nM for REV-ERB $\alpha$ and REV-ERB $\beta$ , respec...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>23</sub> H <sub>31</sub> ClN <sub>4</sub> O <sub>3</sub> S
<b>MW</b>	479.04
<b>CAS Number</b>	[1379686-29-9]
<b>SMILES</b>	<chem>O=C(N1CC(CN(CC2=CC=C(Cl)C=C2)CC3=CC=C([N+])([O-])=O)S3)CC1)NCCCCC</chem>
<b>Expiration Date</b>	12 months from date of receipt.