

Product Datasheet

QCC-374 (orb422383)

Description	QCC-374 is prostanoid agonist potentially for the treatment of pulmonary arterial hypertension.
Note	For research use only
Formula	$C_{28}H_{33}N_3O_2$
MW	443.59
CAS Number	[1356331-63-9]
SMILES	<chem>O=C(O)CCCCCN1CCCC2=NC(C3=CC=C(C)C=C3)=C(C4=CC=C(C)C=C4)N=C21</chem>
Expiration Date	12 months from date of receipt.