



## Product Datasheet

### MEX3B antibody (orb422920)


<b>Description</b>	MEX3B antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC
<b>Immunogen</b>	Recombinant Human RNA-binding protein MEX3B protein (253-375AA)
<b>Target</b>	MEX3B
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>Q6ZN04</b>
<b>Dilution Range</b>	IHC-P:1:200:500, IF/ICC:1:50:200
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of SH-SY5Y C...



Immunohistochemical staining of paraffin...



IHC-P image of human brain tissue using ...