


Product Datasheet

SCUBE1 antibody (orb422935)

Description	SCUBE1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Signal peptide, CUB and EGF-like domain-containing protein 1 protein (453-566AA)
Target	SCUBE1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IWY4
Dilution Range	IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...