



Product Datasheet

SLC16A11 antibody (orb423117)

Description	SLC16A11 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Monocarboxylate transporter 11 protein (428-471AA)
Target	SLC16A11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NCK7
Dilution Range	WB:1:500:5000, IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of A549 cell...



Western blot analysis of A549 whole cel...