

Product Datasheet

Goserelin peptide (orb423181)

Description	Goserelin peptide
Preservatives	The Goserelin peptide was lyophilized with no additives.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Stability: Lyophilized Goserelin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Goserelin should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Please prevent freeze-thaw cycles
Note	For research use only
Application notes	Hormones
Purity	Greater than 98.0% as determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
Solubility (25°C)	It is recommended to reconstitute the lyophilized Goserelin in sterile 18MΩ-cm H ₂ O not less than 100 µg/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	12 months from date of receipt.