



Product Datasheet

IL1 alpha Receptor antibody (orb4233)

Description	IL1 alpha Receptor antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human IL-1R1 (481-569/569 aa)
Target	IL-1R1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	63 kDa
Uniprot ID	P14778
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.




Immunohistochemical staining of rat coli...



Flow cytometric analysis of spleen cell

...



Immunohistochemical staining of rat sple...