

Product Datasheet

Dengue Control antibody (orb423400)

| | |
|--------------------------|---|
| Description | Mouse monoclonal antibody to Dengue Control |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Preservatives | 10mM phosphate buffered saline pH-7.2. |
| Form/Appearance | Sterile Filtered clear colorless solution. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Applications: Immunoassay |
| Clonality | Monoclonal |
| Purity | 90% pure (12% SDS-PAGE and Coomassie staining). |
| Expiration Date | 12 months from date of receipt. |