

Product Datasheet

HA-tag antibody (orb423497)

Description	Mouse monoclonal antibody to HA-tag
Conjugation	Unconjugated
Immunogen	HA-tag conjugated proteins
Preservatives	1mg/ml in PBS (after reconstitution).
Storage	Storage Procedures: In lyophilized form, for long periods, store at 4°C in a dry environment. After reconstitution, if not intended for use within a month, aliquot and store at -20°C
Note	For research use only
Application notes	Titer: By direct ELISA, 1 : 20,000 dilution will yield > 1 O.D using alkaline phosphatase conjugated rabbit anti-mouse Ig (Jackson Laboratories)
Isotype	Mouse IgG2b/kappa
Clonality	Monoclonal
Solubility (25°C)	Reconstitute with H2O to get a 1 mg/ml concentration. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.
Expiration Date	12 months from date of receipt.