

## Product Datasheet

### HBsAg Paired antibody (orb423501)

<b>Description</b>	Mouse monoclonal antibody to HBsAg Paired
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	* HBsAg Labeling antibody: in PBS, pH 7.4.* HBsAg Coating antibody: iin PBS, pH 7.4.
<b>Form/Appearance</b>	2 vials of sterile Filtered clear colorless solution.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Monoclonal Antibodies, Anti-Viral
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Greater than 90%.
<b>Expiration Date</b>	12 months from date of receipt.