

Product Datasheet

Hepatitis B Core antibody (orb423512)

| | |
|--------------------------|---|
| Description | Mouse monoclonal antibody to Hepatitis B Core |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Immunogen | r.HBcAg Core ayw Antigen |
| Preservatives | The protein contains PBS. |
| Form/Appearance | Sterile Filtered solution. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Applications: ELISA, WB, IHC |
| Isotype | Mouse IgG2b |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |