

## Product Datasheet

### HIV-1 gp120 antibody (orb423517)

<b>Description</b>	Rabbit polyclonal antibody to HIV-1 gp120
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Recombinant HIV-1 gp120 produced in baculovirus expression system
<b>Preservatives</b>	Polyclonal Antibody serum was lyophilized in vials.
<b>Form/Appearance</b>	Sterile filtered pinkish powder.
<b>Storage</b>	Stability: Stable at -20°C. Avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Applications: Binds to recombinant HIV-1 gp120 ELISA and Western Blot Optimal dilutions should be determined by the customer, based on titration. Suggested starting dilutions are listed below: Western Blot: 1-500 to 1-1000. ELISA: 1:500 to 1000
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.