

Product Datasheet

Influenza A H5N1 antibody (orb423595)

| | |
|--------------------------|---|
| Description | Rabbit polyclonal antibody to Influenza A H5N1 |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Immunogen | Recombinant Influenza A/H5N1 |
| Preservatives | PBS pH-7.0, 0.005% NaN ₃ and 50% glycerol. |
| Form/Appearance | Sterile Filtered clear solution. |
| Storage | Storage Procedures: Store at -20°C, one month at 4°C |
| Note | For research use only |
| Application notes | Applications: Influenza A haemagglutinin H1N1 Western Blotting at a concentration of 0.5µg/ml. ELISA to be determined |
| Clonality | Polyclonal |
| Purity | 90%. |
| Expiration Date | 12 months from date of receipt. |