

Product Datasheet

Influenza B antibody (orb423596)

Description	Rabbit polyclonal antibody to Influenza B
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	Recombinant Influenza B Jiangsu/10/2003
Preservatives	PBS pH 7.0, 0.005% Sodium Azide & 50% glycerol.
Storage	Storage Procedures: Store at -20°C
Note	For research use only
Application notes	Applications: Influenza B haemagglutinin Western Blotting at a concentration of 0.5µg/ml. ELISA to be determined
Clonality	Polyclonal
Purity	90%.
Expiration Date	12 months from date of receipt.