

## Product Datasheet

### Influenza-A H1N1 antibody (orb423597)

<b>Description</b>	Mouse monoclonal antibody to Influenza-A H1N1
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Influenza A hemagglutinin H1N1
<b>Preservatives</b>	In PBS, pH 7.4 and 0.09% NaN <sub>3</sub> .
<b>Form/Appearance</b>	Sterile Filtered clear solution.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Applications: Influenza A haemagglutinin 1 H1N1 immunodetection in direct or indirect ELISA, sandwich immunoassay, Western Blotting. Recommended pairs for Influenza A H1N1 in sandwich immunoassay (coating - conjugate): All pairs detect H1N1 of Influenza A as well as recombinant H1N1 with high specificity and do not interact with haemagglutinin H3N2
<b>Isotype</b>	mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.