

## Product Datasheet

### Influenza-A Paired antibody (orb423602)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Mouse monoclonal antibody to Influenza-A Paired   |
| <b>Reactivity</b>        | Virus   |
| <b>Conjugation</b>       | Unconjugated  |
| <b>Preservatives</b>     | * Influenza-A gold conjugation antibody (7.4mg/1ml) in PBS, pH 7.4. * Influenza-A capture antibody (6.4mg/1ml) in PBS, pH 7.4.          |
| <b>Form/Appearance</b>   | 2 vials of sterile Filtered clear colorless solution.   |
| <b>Storage</b>           | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Applications: Lateral flow immunoassay  |
| <b>Clonality</b>         | Monoclonal  |
| <b>Purity</b>            | Greater than 90%.   |
| <b>Expiration Date</b>   | 12 months from date of receipt.   |