

Product Datasheet

Influenza-B Paired antibody (orb423604)

Description	Mouse monoclonal antibody to Influenza-B Paired
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	* Influenza-B conjugation antibody in PBS, pH 7.4 . * Influenza-B coating antibody in PBS, pH 7.4 .
Form/Appearance	2 vials of sterile Filtered clear colorless solution.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: Lateral flow immunoassay
Clonality	Monoclonal
Purity	Greater than 90%.
Expiration Date	12 months from date of receipt.