

Product Datasheet

M. Tuberculosis RV1681 antibody (orb423650)

Description	Rabbit polyclonal antibody to M. Tuberculosis RV1681
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	Glycine-Tris buffer, pH 7.5.
Form/Appearance	Sterile Filtered clear colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Polyclonal Antibodies, Anti-Viral
Clonality	Polyclonal
Purity	Pure Rabbit IgG, > 95% pure as determined by 12% PAGE (coomassie staining).
Expiration Date	12 months from date of receipt.