

Product Datasheet

Norovirus Group-I Paired antibody (orb423708)

Description	Mouse monoclonal antibody to Norovirus Group-I Paired
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	* Norovirus Group I capture antibody in 1xPBS, pH 7.4.* Norovirus Group I conjugating antibody in 1xPBS, pH 7.4.
Form/Appearance	2 vials of sterile filtered clear colorless solution.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: Lateral flow immunoassay
Clonality	Monoclonal
Purity	Greater than 95%.
Expiration Date	12 months from date of receipt.