

Product Datasheet

Norovirus Group-II Paired antibody (orb423709)

| | |
|--------------------------|---|
| Description | Mouse monoclonal antibody to Norovirus Group-II Paired |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Preservatives | * Norovirus Group II capture antibody in 1xPBS, pH 7.4.* Norovirus Group II conjugating antibody in 1xPBS, pH 7.4. |
| Form/Appearance | 2 vials of sterile filtered clear colorless solution. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Applications: Lateral flow immunoassay |
| Clonality | Monoclonal |
| Purity | Greater than 90%. |
| Expiration Date | 12 months from date of receipt. |