

## Product Datasheet

### S.typhi Paired antibody (orb423805)

<b>Description</b>	Mouse monoclonal antibody to S.typhi Paired
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	S.Typhi capture&conjugate antibody (2.7mg/1ml) in 0.01M phosphate buffered saline, pH 7.2 and 0.1% sodium azide.
<b>Form/Appearance</b>	Sterile filtered clear colorless solution.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Monoclonal Antibodies, Anti-Viral
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Greater than 90%.
<b>Expiration Date</b>	12 months from date of receipt.