

Product Datasheet

Shiga Like Toxin 1B antibody (orb423815)

Description	Mouse monoclonal antibody to Shiga Like Toxin 1B
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	1x PBS and 0.05% sodium nitrate.
Form/Appearance	Sterile Filtered clear colorless solution.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Monoclonal Antibodies, Anti-Viral
Clonality	Monoclonal
Purity	Greater than 95% as determined by 12% PAGE (coomassie staining).
Expiration Date	12 months from date of receipt.