

Product Datasheet

Zika NS1 Paired antibody (orb423895)

Description	Mouse monoclonal antibody to Zika NS1 Paired
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	* Zika NS1 gold conjugation antibody in PBS, 10mg rAlbumin/ml & 0.1% Sodium Nitrate.* Zika NS1 capture antibody in PBS, 10mg rAlbumin/ml & 0.1% Sodium Nitrate.
Form/Appearance	2 vials of sterile filtered clear colorless solution.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: Lateral flow immunoassay
Clonality	Monoclonal
Purity	Greater than 95%.
Expiration Date	12 months from date of receipt.