

Product Datasheet

Human ADPRH Protein (orb423994)

Description	Recombinant of human ADPRH protein
Preservatives	The ADPRH solution (0.5mg/ml) contains 20mM Tris-HCl buffer, pH 8.0, 10% glycerol, 1mM DTT and 100mM NaCl.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SGLVPRGSH MGSHEKYVA AMVLSAAGDA LGYNGKWEF LQDGEKIHQR LAQLGGLDAL DVGRWRVSDD TVMHLATAEA LVEAGKAPKL TQLYLLAKH YQDCMEDMDG RAPGGASVHN AMQLKPGKPN GWRIPFNSHE GGCGAAMRAM CIGLRFPHHS QLDTLIQVSI ESGRMTHHHP TGYLALASA LFTAYAVNSR PPLQWGKGLM ELLPEAKKYI VQSGYFVEEN LQHWSYFQTK WENYLKLRGI LDGESAPTFP ESFGVKERDQ FYTSLSYSGW GGSSGHDAPM IAYDAVLAAG DSWKELAHRA FFHGGSDSDST AAIAGCWWGV MYGFKGVSPS NYEKLEYRNRLEETARALYS LGSKEDTVIS L
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
Expiration Date	12 months from date of receipt.