

Product Datasheet

Human APOBEC4 Protein (orb424126)

Description	Recombinant of human APOBEC4 protein
Preservatives	APOBEC4 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.4M UREA and 10% glycerol.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please avoid freeze thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMEPIYEE YLANHGTVK PYYWLSFSLD CSNCPYHIRT GEEARVSLTE FCQIFGFPYG TTFPQTKHLT FYELKTSSGS LVQKGHASSC TGNVIHPESM LFEMNGYLDL AIYNNDSIRH IILYSNNSPC NEANHCCISK MYNFLTYPG ITLSIYFSQL YHTEMDFPAS AWNREALRSL ASLWPRVVLV PISGGIWHVSV LHSFISGVSG SHVFQPILTG RALADRHAY EINAITGVKP YFTDVLLQTK RNPNTKAQEA LESYPLNNAF PGQFFQMPSG QLQPNLPPDL RAPVVFVLP LRDLPMMHMG QNPKNPRNIV RHLNMPQMSF QETKDLGRLP TGRSVEIVEI TEQFASKEA DEKKKKKGGK
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.