

Product Datasheet

Human ATOX1 Protein (orb424209)

Description	Recombinant of human ATOX1 protein
Preservatives	ATOX1 solution containing 20mM Tris pH-8, 1mM DTT and 10% glycerol.
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MPKHEFSVDM TCGGCAEAVS RVLNKLGGVK YDIDLPNKKV CIESEHSMDT LLATLKKTKG TVSYLGLE
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.